

A-Scan


Touch Screen Table Top

OPTISCAN



Diagnosis	Features	Aesthetics
<ul style="list-style-type: none"> • Axial Length measurement • Depth of anterior chamber • Thickness of lens • Standard deviation 	<ul style="list-style-type: none"> • Calculating IOL power using 6 formulae • Built-in thermal printer • 100 patients data • Auto delete wrong readings with Artificial Intelligence • USB port for PC printers 	<ul style="list-style-type: none"> • Compact and light weight • Ultra light weight probe • Easy to use touch screen

A-SCAN SPECIFICATION



Clinical Accuracy:

- New enhanced algorithms for rapid measurement
- 10 measurements presented on a single screen
- Multiple velocities and lens type for higher accuracy
- Statistical analysis of axial length, vitreous, anterior chamber depth and lens thickness
- Continuous and manual readings

Data Analysis:

- All current generation IOL power calculation formulae
- View multiple power calculations using different formulae in a single screen

Transducer:

- Frequency : 10MHz +/- 10%
- Clinical accuracy : +/- 0.1 mm
- Theoretical accuracy : +/- 0.05 mm
- Axial length : 15 - 40 mm
- Gain range : 100 dB
- Fixation light : Internal fixation LED

IOL Power Calculation Formulae:

- | | |
|-------------|-----------------|
| 1) SRK T | 4) Binkhorst II |
| 2) SRK II | 5) Hoffer-Q |
| 3) Holladay | 6) Haigis STD |

Eye Modes:

- | | |
|------------|-------------|
| 1) Normal | 5) Silicone |
| 2) Dense | 6) Acrylic |
| 3) Aphakic | 7) Custom |
| 4) PMMA | |

Standard Accessories:

- Solid tip transducer with fixation light
- Built-in test cylinder
- User manual
- Foot pedal
- Printer paper - 4 Rolls
- Power cable

Other:

- Measurement mode : Contact and Immersion
- Print out : Thermal print out
- Power supply : Universal
- Data feeding / Menu selection : 7" touch screen

Salient Features




OPTITECH[®]
 TARUN ENTERPRISES
 An EN ISO 13485:2016 Company

H.O.: 8/8, Strachy Road, Prayagraj (Allahabad) – 211001, U.P., INDIA
 Phone: +91 - 8176080204, Email: info@optitecheyecare.com

B.O.: A96/1, Second Floor, Phase II, Mayapuri Industrial Area, New Delhi – 110064
 Phone: 011-47070362, Email: delhi@optitecheyecare.com

www.optitecheyecare.com +91 - 9335154556 OptitechOfficial

Your Nearest Dealer



*Specifications are subject to change without notice